
## **thermo**scientific





Thermo Scientific™ Benchtop Shakers are robust enough for everyday use yet provide precise sample agitation for all your applications. With a broad range of digital speed and time settings, you can be confident that your sample will shake, mix, or pulse at the desired rate every time.

Regardless of your application need, you'll be able to set your sample agitation parameters and move on to your next task, helping to bring increased efficiency to your lab.

### Product highlights

- Models designed for a variety of controlled agitations, including shaking, rocking, waving rotation, touch vortex mixing, and more.
- Newer models come standard with digital controls and upgrades that allow for increased sample mixing
- Adjustable tilt angles by digital control no additional tools needed
- Robust body and standard rubber feet ensure shakers stay in place on the lab bench





# Contents

Shaker	
Digital Microplate Shaker	4
Mixers	
Digital Vortex Mixer	6
MaxiMix Vortex Mixer	8
Rotators	
Digital Waving Rotator	10
Digital Mini Rotator	12
MultiTube Rotators	14
Cel-Gro Tissue Culture Rotator	16
Rockers	
Digital Rocker	18
Vari-Mix Platform Rocker	20
Vari-Mix and Speci-Mix Test Tube Rockers	22
Rollers	
Digital Bottle/Tube Roller	24
Digital Tube Revolver	26



Thermo Scientific™ Digital Microplate Shakers are developed for precise and accurate shaking of microplates for immunoassay, screening, magnetic bead, and other applications.

### **Key Features**

• Orbit diameter: 2.5mm

• Digital speed and time displays

• Speed range: 150 to 1,200rpm

• Speed control:

- ≤299rpm: ±1rpm

- 300 to 1,200rpm: ±2rpm

Warranty: 2 years | Certifications: cCSAus, CE

• Timer range: up to 99 hr. 59 min

• Large platform capacity: 6 x 96-well plates

• Maximum load: 4kg (lb.)




Digital Microplate Shaker shown with Thermo Scientific Clamp for 96-well Microplates.

Product	Electrical Requirements	Tray Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Microplate Shaker	100 to 240V. 50/60Hz	27.5 x 23.5cm	38.3 x 26 x 14.5cm (15.1 x 10.2 x 5.7 in.)	12.8kg	120V USUP	88882005
Digital Microplate Shaker	100 to 240V, 50/60HZ	(10.8 x 9.3 in.)		(28.2 lb.)	230V EU/UK/CHN	88882006

### **Accessories**

Product and Description	Quantity	Catalog Number
Large Platform (holds 6 microplates)	Each	88882114
Clamp for 96-well microplate	Each	88882115
Large Assembly Platform	Each	88882116
Tube Rack, 50 x 11mm (1.5mL)	Each	88882117
Tube Rack, 40 x 18mm (1.5mL)	Each	88882118
Tube Rack, 21 x 30mm (50mL)	Each	88882119
Screw, GB/T819.1, 4 x 8, 4/pack	Each	88882130
Screw, GB/T819.1, 5 x 10, 4/pack	Each	88882131
Screw, GB/T9074.4, 3 x 6, 2/pack	Each	88882132
Power Adapter with US plug	Each	88870126
Power Adapter with AUS/CN plug	Each	88870127
Power Adapter with EU plug	Each	88870128
Power Adapter with UK plug	Each	88870129

Includes standard platform with 6 clamps (88882115) installed; holds 4 microplates.



Digital Microplate Shaker accessories shown include Large Assembly Platform and Large Platform that hold up to 6 microplates (top), and Tube Racks holding 21, 40, and 50 tubes respectively (bottom).



The Thermo Scientific™ Digital Vortex Mixer is designed for a wide variety of applications.

### **Key Features**

- Digital speed and time displays
- Continuous and touch modes
- Speed range: 200 to 3,000rpm
- Speed accuracy:
  - 200 to 1,500rpm: ±5%
  - 1,500 to 3,000rpm: ±10%

- Rapid acceleration to maximum rpm
  - Touch: 1,500rpm in 3 seconds
  - Continuous: 3,000rpm in 6 seconds
- Countdown timer (up to 99 minutes)
- Maximum load: 0.5kg at 800rpm



Includes rubber platform (Cat. No. 88882120). Other accessories sold separately. Warranty: 2 years | Certifications: cCSAus, CE

Product	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Vortex Mixer		20.5 x 15.1 x 13.1cm		120V USUP	88882009
Digital Vortex Mixer	100 to 240V, 50/60Hz	(8.1 x 5.9 x 5.2 in.)	6.1kg (13.5 lb.)	230V EU/UK/ CHN	88882010

Digital Vortex Mixer shown with platform for flat-bottom vessels.

Product and Description	Mode	Recommended Maximum Speed	Quantity	Catalog Number
Rubber Platform	Touch/continuous	3,000rpm	Each	88882120
Platform for flat-bottomed vessels	Touch/continuous	3,000rpm	Each	88882121
Microplate Tray	Continuous	1,200rpm	Each	88882122
Tube Holder: 19 x 10mm	Continuous	2,000rpm	Each	88882123
Tube Holder: 13 x 12mm	Continuous	1,800rpm	Each	88882124
Tube Holder: 12 x 15mm	Continuous	1,000rpm	Each	88882125
Tube Holder: 7 x 20mm	Continuous	1,000rpm	Each	88882126
Tube Holder: 4 x 26mm	Continuous	900rpm	Each	88882127
Multi-Power Adapter, 220V	_	_	Each	88870124
Power Adapter, 120V	_	_	Each	88870125



Vortex Mixer Plugs and Optional Platforms



The Thermo Scientific™ MaxiMix™ Vortex Mixer is designed for enzymatic and RIA assays, general test tube mixing, biochemical assays, viral dilutions preparations, precipitated assays, cell suspension vortexing, and tissue sample mixing.

### **Key Features**

- For single or multiple tubes or small closed flasks
- Simultaneously mixes up to four test vessels
- Top-mounted push-button ON/OFF switch
- Vary vortex speed by changing the tube pressure against the foam rubber top
- Durable, white plastic case resists acids and alkalis
- Cast-aluminum base with suction-cup rubber feet for stability
- Compact design; 4 in. diameter (10.2cm) foam rubber top
- Speed up to 3,000rpm



MaxiMix Vortex Mixer

Includes one extra foam top (PT167X2), 3-wire power cord, and plug.

Warranty: 90 days labor; one year parts and labor | Certifications: CE listed (230V models)

Product	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
MaxiMix Vortex Mixer	120-240V 50/60Hz	10.0 v 10.0 v 10.6cm (F.F. v F.F. 4.0 in.)	1 Elea (2 E lb.)	120V USUP	M16715Q
MaxiMix Vortex Mixer	120-240V 50/60Hz	13.9 x 13.9 x 10.6cm (5.5 x 5.5 4.2 in.)	1.5kg (3.5 lb.)	230V EU/UK/CHN	M16710-33Q



The Thermo Scientific<sup>™</sup> Digital Waving Rotator is designed for ELISA, protein synthesis, gel staining and de-staining gels, antibody staining, and other time-sensitive applications.

### **Key Features**

• Digital speed and time displays

• Adjustable angle: 1° to 12°

• Speed range: 5 to 120rpm

• Speed control:

≤30rpm: ±1rpm>30rpm: ±2rpm

• Timer range: up to 99 hr. 59 min

• Capacity: 6 x 96-well plates

• Maximum load: 9.8kg (22 lb.)



Includes Dimpled Rubber Mat (88882102), Rubber Strips (88882103, 6/Pack). Use optional Rubber Mat (88882101) for plates, dishes, culture bottles, and bags. Warranty: 2 years | Certifications: CE

Digital Waving Rotator shown with Thermo Scientific™ Nunc™ Test Tubes on a Thermo Scientific Rubber Mat. Nunc Test Tubes sold separately.

Product	Electrical Requirements	Tray Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digtial Waving Rotator	100 to 240V, 50/60Hz	31.6 x 27.6 x 1.9cm	38.3 x 29.5 x 22.4cm	13.5kg	120V USUP	88882003
Digtial Waving Rotator		(12.4 x 10.9 x 0.7 in.)	(15.1 x 11.6 x 8.8 in.)	(29.8 lb.)	230V EU/UK/CHN	88882004

Product and Description	Quantity	Catalog Number
Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	88882101
Dimpled Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	88882102
Rubber Strips	6/pack	88882103
Replacement Rubber Tube	Each	88882128
Screw for Rubber Strips	8/Pack	88882129
Power Adapter with US plug	Each	88870126
Power Adapter with AUS/CN plug	Each	88870127
Power Adapter with EU plug	Each	88870128
Power Adapter with UK plug	Each	88870129



Waving Rotator accessories shown include: Replacement Rubber Tube and Rubber Strips (upper left), Screw for Rubber Strips (upper right), Rubber Mat and Dimpled Rubber Mat (lower left and right).



The Thermo Scientific™ Digital Mini Rotator is designed for culturing bacteria and fungi, hybridization assays, and gel staining and de-staining.

### **Key Features**

Orbit diameter: 10mm

• Speed range: 50 to 500rpm

• Speed control:

- ≤100rpm: ±1rpm - 101 to 500rpm: ±1%

- Continuous operation or timer modes
- Maximum load: 4.5kg at 300rpm



Includes Universal Aluminum Platform (88882105). Multi-platform options available for loading vessels of all sizes. Warranty: 2 years | Certifications: cCSAus, CE

Digital Mini Rotator shown with Thermo Scientific Flask Clamp – set of 6 500mL clamps with flasks. Flasks sold separately.


Product	Electrical Requirements	Tray Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Mini Rotator	100 to 040 / 50/00  -	31 x 28.4cm	38.3 x 28.4 x 13cm	14.4kg	120V USUP	88882007
Digital Mini Rotator	100 to 240V, 50/60Hz	(12.2 x 11.2 in.)	(15.1 x 11.2 x 5.1 in.)	(31.8 lb.)	230V EU/UK/CHN	88882008


### **Accessories**

Product and Description	Quantity	Catalog Number
Universal Attachment	Each	88882104
Universal Platform, aluminum	Each	88882105
Flack Clamp, 50mL, max: 20 per platform	Each	88882106
Flask Clamp, 100/125mL, max: 8 per platform	Each	88882107
Flask Clamp, 250mL, max: 8 per platform	Each	88882108
Flask Clamp, 500mL, max: 6 per platform	Each	88882109
Flask Clamp, 1000mL, max: 4 per platform	Each	88882110
Flask Clamp, 2000mL, max: 2 per platform	Each	88882111
Sticky Mat, 14 x 14cm	Each	88882112
Flat Platform with rubber mat	Each	88882113
Tube Rack, 50 x 11mm (1.5mL tubes)	Each	88882117
Tube Rack, 40 x 18mm (1.5mL tubes)	Each	88882118


Product and Description	Quantity	Catalog Number
Tube Rack, 21 x 30mm (50mL tubes)	Each	88882119
Screw, GB/T819.1, 4 x 8	4/pack	88882130
Screw, GB/T819.1, 5 x 10	4/pack	88882131
Screw, GB/T9074.4, 5 x 10	6/pack	88882133
Power Adapter with US plug	Each	88870126
Power Adapter with AUS/CN plug	Each	88870127
Power Adapter with EU plug	Each	88870128
Power Adapter with UK plug	Each	88870129







Digital Mini Rotator accessories shown include (from left to right) Universal Attachment, Clamps, Universal Platform.



Thermo Scientific™ MultiTube Rotators are designed for hematology, blood chemistry, DNA extractions, and other laboratory applications where a continuous, gentle motion is necessary.

### **Key Features**

- Circular motion and combined rocking pattern reduces the potential for shearing or other damage
- Fixed speed: 30rpm
- Rotation angle: variable 20° from horizontal
- Compact size can be used in smaller spaces
- Portable, lightweight metal body is durable and easy to clean


- Four rubber feet for added stability
- No need for racks or clamps platform side walls and removable dimpled silicone pad (included) hold tubes in place
- Temperature and humidity range: ambient: 15° to 40°C (59° to 104°F); 20 to 80% RH non-condensing



MultiTube Rotator shown with test tubes (sold separately).

Warranty: 90 days labor (North America only); one year parts | Certifications: 240V model is CE listed

Product	Tube Capacity	Load Capacity	Electrical Requirements	Mixing Surface	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
MultiTube Rotator	14 x 7. 10	1.4 kg	120V 60Hz, 0.3A	13 x 22cm	21 x 22 x 18cm	3.6kg	120V USUP	4632Q
MultiTube Rotator	or 15mL			(8.5 x 8.8 x 8.1 in.)	(8 lb.)	230V EU/UK/ CHN	4632-1CEQ	



The Thermo Scientific™ Cel-Gro Tissue Culture Rotator is designed for cell culture.

### **Key Features**

- Adjustable angle and speed
- Durable and corrosion-resistant powder-coated stainless-steel outer body construction
- Adjust elevation angle using a screw on the base

- Pilot light indicates when unit is ON
- Speed range: 12 to 70rpm
- Temperature and humidity range: up to 40°C (104°F); 20 to 80% non-condensing



Warranty: 90 days labor (North America only); one year parts | Certifications: cCSAus approved

Product	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Cel-Gro Tissue Culture Rotator	120V 60Hz, 0.5A	38.1 x 35.6 x 57.2cm (15 x 14 x 22.5 in.)	9kg (19 lb.)	120V USUP	1640Q

Cel-GroTissue Culture Rotator assembled with Thermo Scientific Rotator Drum Accessory and tubes (sold separately).

Product	Capacity	Catalog Number
Rotator Drum	142 x 17mm tubes	1647Q
Rotator Drum	76 x 26mm tubes	1648Q
Rotator Drum	60 x 30mm tubes	1651Q



The Thermo Scientific<sup>™</sup> Digital Rocker is designed for ELISA, protein synthesis, gel staining and de-staining, antibody staining, and other time-sensitive applications.

### **Key Features**

• Digital speed and time displays

• Adjustable rocking angle: 1° to 15°

• Speed range: 1 to 70rpm

Speed control:
- ≤50rpm: ±1rpm
- >50rpm: ±2rpm

• Timer range: continuous or 1 min. to 99 hr. 59 min.

• Capacity: 6 x 96-well plates

• Maximum load: 4.5kg at ≤50rpm




Includes Dimpled Rubber Mat (88882102) and Rubber Strips (88882103, 6/Pack). Optional non-slip rubber mat tray is designed for carrying plates, dishes, culture bottles, and bags. Warranty: 2 years | Certifications: cCSAus, CE

Digital Rocker shown with flasks on top of a Thermo Scientific Dimpled Rubber Mat. Flasks sold separately.

Product	Electrical Requirements	Tray Dimensions D x W x H	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Rocker	100 to 240V. 50/60Hz	32.2 x 28.2 x 1.9cm	39 x 29.2 x 20.6cm	13.5kg	120V USUP	88882001
Digital Rocker	100 to 240V, 50/60M2	(12.7 x 11.1 x 0.7 in.)	(15.4 x 11.5 x 8.1 in.)	(29.8 lb.)	230V EU/UK/CHN	88882002

Product and Description	Quantity	Catalog Number
Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	88882101
Dimpled Rubber Mat, 32.2 x 28.2cm (12.7 x 11.1 in)	Each	88882102
Rubber Strips	6/pack	88882103
Power Adapter with US plug	Each	88870126
Power Adapter with AUS/CN plug	Each	88870127
Power Adapter with EU plug	Each	88870128
Power Adapter with UK plug	Each	88870129



Thermo Scientific Rubber Strips (top), Rubber Mat (bottom left), and Dimpled Rubber Mat (bottom right).



The Thermo Scientific<sup>™</sup> Vari-Mix<sup>™</sup> Platform Rocker is designed for applications such as hybridization, blotting and staining, or destaining gels.

### **Key Features**

- Built-in timer
- Large, white nonskid rubber mixing platform allows easy content viewing and removes easily for cleaning
- Easy-to-clean chemical-resistant finishes
- Optional double-tier platform for increased capacity
- Speed range: 5 to 30rpm

- Rocking angle: 1° to 48°
- Timer range: up to 2 hr. (60Hz); up to 3 hr. (50Hz)
- Load capacity: 6.8kg (15 lb.)
- Temperature range: 4° to 40°C (39.2° to 104°F)
- Operation humidity range: 20 to 80% non-condensing



Vari-Mix Platform Rocker with Double-Tier Platform attached (flasks and containers sold separately).

Includes 3-wire cord and plug; 240V	model also includes European cord set.
Warranty: 90 days labor (North Ame	rica only); one year parts

Product	Electrical Requirements	Platform Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Vari-Mix Platform Rocker	120V 50/60Hz, 0.10A, 12W	32 x 26.9cm (12.6 x 10.6 in.)	30.4 x 39.6 x 18.5cm (12 x 15.6 x 7.3 in.)	2.7kg (6 lb.)	US	M79735Q

Product	Description	Catalog Number
Double-Tier Platform	32 x 26.9cm (12.6 x 10.6 in.)	AY797X1



Thermo Scientific™ Vari-Mix™ and Speci-Mix™ Test Tube Rockers are designed for uniform specimen suspension.

### **Key Features**

- Choose fixed or adjustable speeds
- Hold tubes with diameters of 10 to 30mm, including blood collection tubes
- Reversible lipped white silicon pad (included); one side for 110mm (4.33 in.) tubes, longer tubes on reverse
- Autoclavable white silicone rubber pad for enhanced viewing
- Temperature and humidity range: ambient: 4° to 40°C (39.2° to 104°F); 20 to 80% RH non-condensing



Warranty: 90 days labor (North America only); one year parts | Certifications: CE listed (240V models only)

Product	Capacity	Speed	Rocking Angle	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Speci-Mix Test	8 x 10 to 20mm;	18rpm,	48°, fixed	120-240V 50/60Hz	13.9 x 19.1 x 10.7cm	1 (1/0 /0 lb.)	120V USUP	M71015Q
Tube Rocker	3 x 22 to 30mm	fixed	46 , lixea	.20 2 10 1 00, 001 .2	(5.5 x 7.5 x 4.3 in.)	1.4kg (3 lb.)	230V EU/UK/CHN	M71010-33Q
Vari-Mix Test	16 x 10 to 20mm;	2 to	Adjustable: 1° to 48°	120-240V 50/60Hz	13.9 x 38.1 x 10.7cm	2.6kg (5.8 lb.)	120V USUP	M48725Q
Tube Rocker	7 x 22 to 30mm	20rpm	Adjustable. 1 to 46	120 2 10 7 00, 001 12	(5.5 x 15 x 4.3 in.)	2.0kg (5.0 lb.)	230V EU/UK/CHN	M48720-33Q



The Thermo Scientific<sup>™</sup> Digital Bottle/Tube Roller is designed for cell culture, hybridization assays, liquid mixing, and bioproduction.

### **Key Features**

- LED roller speed display
- Roller speed: 1 to 80rpm
- Number of rollers: six standard; add up to 11
- Timing range: continuous or 0 to 999 minutes
- Optional remote control operation
- Anti-corrosive and wear-resistant construction for use in humid and CO<sub>2</sub> environments
- Stacking bars let you stack two units high to save valuable bench space



Includes O-rings for smooth rolling of 15 or 50mL tubes. Warranty: 2 years | Certifications: CE, cCSAus

Two Digital Bottle/Tube Rollers assembled with accessory Bars for Stacking. Thermo Scientific Nunc Tubes and Bottles (displayed) sold separately.

Product	Electrical Requirements	Roller Dimensions	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Bottle/Tube Roller	100V-240V. 50Hz 60Hz	28 x 259mm	378 x 360 x 72mm	6 Elva (14 0 lb.)	US	88881003
Digital Bottle/Tube Roller	1000-2400, 5002 6002	(1.1 x 10.2 in.)	(14.9 x 14.2 x 2.8 in.)	6.5kg (14.3 lb.)	Various	88881004

Product and Description	Quantity	Catalog Number
Additional Rollers	2/pack	88881010
O-Ring for 50mL Tubes, 23.5 x 3.55mm (0.93 x 0.14 in.)	20/pack	88881011
O-Ring for 15mL Tubes, 13.2 x 2.65mm (0.52 x .1 in.)	20/pack	88881012
Bars for Stacking	4/pack	88881013
Remote Control	Each	88881014





The Thermo Scientific™ Digital Tube Revolver is designed for molecular biology, clinical, and biochemistry applications.

### **Key Features**

- Interchangeable rotisseries to accommodate different tube sizes
- Paddles can be adjusted from 0 to 90 degrees
- Choose agitation mode; agitation will occur at every 40 degrees in the rotation cycle
- Small footprint, plug-and-play simplicity
- Maintenance-free motor, quiet operation
- Rotation speed: 10 to 40rpm
- Operating temperature: 4 to 60°C



Digital Tube Revolver shown with Paddle for 10 or 15mL tubes (front) and Paddle for 21 0.5mL or 0.8mL tubes (rear). Thermo Scientific Nunc Tubes (displayed) sold separately.

### Paddles included:

- For 7 x 10/15mL tubes (2)
- For 7 x 5mL/7mL tubes (2)
- For 15 x 1.5mL/2mL tubes (2)
- For 21 0.5mL/0.8mL tubes (2)
- For 6 x 50mL tubes (1)

Warranty: 2 years | Certifications: cCSAus, CE

Product	Electrical Requirements	Overall Dimensions D x W x H	Weight	Plug Type	Catalog Number
Digital Tube Revolver, 120V	100/ 040/ 50/1- 00/1-	225 x 154 x 195mm	1.01.01 (0.05 lb.)	US	88881001
Digital Tube Revolver, 120V	100V-240V, 50Hz-60Hz	(8.9 x 6 x 7.7 in.)	1.2kg (2.65 lb.)	Various	88881002

Product and Description	Quantity	Catalog Number
Paddles for 10mL/15mL tubes	2/pack	88881005
Paddles for 5mL/7mL tubes	2/pack	88881006
Paddles for 1.5mL/2mL tubes	2/pack	88881007
Paddles for 0.5mL/0.8mL tubes	2/pack	88881008
Paddle for 50mL tubes	Each	88881009



Digital Tube Revolver Paddles



## Find out more at **thermofisher.com/benchtopshakers**

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific use or application. © 2012-2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. B-SHAKEROTATEMIX-E 0220

